



# **Progenesis LC-MS User Guide**

**Analysis workflow guidelines**

**for version 4.1**



## Contents

Introduction .....	3
How to use this document .....	3
How can I analyse my own runs using LC-MS? .....	3
LC-MS Data used in this user guide .....	3
Workflow approach to LC-MS run analysis .....	4
Restoring the LC-MS Tutorial .....	5
Stage 1: Import Data and QC review of LC-MS data set .....	6
Stage 2: Automatic Alignment of your data .....	7
Stage 3: Licensing .....	9
Stage 4: Review Alignment .....	9
Stage 5A: Filtering .....	12
Stage 5B: Reviewing Normalisation .....	16
Stage 6: Experiment Design Setup for Analysed Runs .....	19
Stage 7: Review Peak Picking and editing of results .....	22
Stage 8: Peptide Statistics on selected features .....	31
Stage 9: Identify peptides .....	35
Stage 10: Refine Identifications .....	39
Stage 11: Resolve Conflicts .....	40
Stage 12: Review Proteins .....	45
Stage 13: Protein Statistics .....	48
Stage 14: Reporting .....	49
Creating an Inclusion list .....	52
Congratulations! .....	53
Appendix 1: Stage 1 Data Import and QC review of LC-MS data set .....	54
Appendix 2: Stage 1 Data QC review and addition of exclusion areas .....	59
Appendix 3: Licensing runs (Stage 3) .....	60
Appendix 4: Manual assistance of Alignment .....	61
Appendix 5: Within-subject Design .....	66
Appendix 6: Power Analysis (Peptide Stats) .....	68
Appendix 7 Search engine parameters (Stage 9) Mascot .....	69

## Introduction

This user guide takes you through a complete analysis of 6 LC-MS runs with 2 groups (3 replicate runs per group) using the unique Progenesis LC-MS workflow. It starts with LC-MS data file loading then Alignment, followed by Analysis that creates a list of interesting features (peptides) which are explored within Progenesis Stats using multivariate statistical methods then onto Protein identity and Reporting.

To allow ease of use the tutorial is designed to start with the restoration of an Archived experiment where the data files have already been loaded. However, the document covers all the stages in the LC-MS workflow, therefore if you are using your own data files please refer to Appendix 1 (page 54) then start at page 6.

## How to use this document

You can print this user guide to help you work hands-on with the software. The complete user guide takes about 50 minutes (dependant on PC spec) and is divided into two sections. This means you can perform the first half focused on LC-MS run alignment and analysis then complete the second half of analysis exploring comparative differences and Protein identity at a convenient time. If you experience any problems or require assistance, please contact us at [support@nonlinear.com](mailto:support@nonlinear.com)

## How can I analyse my own runs using LC-MS?


You can freely explore the quality of your LC-MS data using Data Import and then licence your own LC-MS runs using this evaluation copy of Progenesis LC-MS. Instructions on how to do this are included in a section at the end of the user guide document. Alternatively if you would like to arrange a demonstration in your own laboratory contact [support@nonlinear.com](mailto:support@nonlinear.com) and we will help you.














## LC-MS Data used in this user guide

NLD would like to thank Dr Robert Parker and Prof Haroun Shah at the Health Protection Agency, London, UK for providing the example data used in this user guide as well as invaluable discussion on the handling of the data.

## Workflow approach to LC-MS run analysis

Progenesis LC-MS adopts an intuitive **Workflow** approach to performing comparative LC-MS data analysis. The following user guide describes the various stages of this workflow (see below) focusing mainly on the stages from Alignment to Report.




Stage	Description	Page
	<b>LC-MS Import Data:</b> Selection and review of data files for analysis.	5
	<b>Automatic Alignment:</b> Automatic Reference selection and alignment	7
	<b>Licensing:</b> allows licensing of individual data files when there is no dongle attached (Appendix 3)	9
	<b>Review Alignment:</b> automatic and manual run alignment	9
	<b>Filtering:</b> defining filters for peaks based on Retention Time, m/z , Charge and Number of Isotopes.	13
	<b>Review Normalisation:</b> explains LC-MS normalisation	16
	<b>Experiment Design Setup:</b> defining one or more group set ups for analysed aligned runs	19
	<b>Review Peak Picking:</b> review and validate results, edit peak detection, tag groups of peaks and select peaks for further analysis	22
	<b>Peptide Statistics:</b> performing multivariate statistical analysis on tagged and selected groups of peptides	31
	<b>Identify Peptides:</b> managing export of MS/MS spectra to, and import of peptide ids from Peptide Search engines	35
	<b>Refine Identifications:</b> manage peptide ids and filters	39
	<b>Resolve Conflicts:</b> validation and resolution of peptide id conflicts for data entered from Database Search engines	40
	<b>Review proteins:</b> review protein and peptide identity	45
	<b>Protein Statistics:</b> multivariate statistical analysis on proteins	48
	<b>Report:</b> generate a report for proteins and/or peptides	49

## Restoring the LC-MS Tutorial


Open Progenesis LC-MS and download the Compressed (.zip) Tutorial Archive file from the the 'View online tutorial' link shown below, placing it in a **new folder** on your desktop. Before restoring the tutorial in the software **you must** first right click on the (.zip) file and extract it to the same folder.

Now you can restore the uncompressed LC-MS tutorial archive file. To do this, first locate the LC-MS Tutorial Archive file using the **Open** button and press Open.




This opens the 'Import from archive' dialog.

Select the **Create a new experiment** option and select the folder in which you placed the archive, using the icon (to the right) .



Then press **Import**.




**Note:** use the **Replace an existing experiment** option if you want to over-write an existing version of the tutorial.

**Tip:** at each stage in the software there are links to more information and help on the website.


## Stage 1: Import Data and QC review of LC-MS data set

The LC-MS tutorial will now open at the LC-MS **Import Data** stage (see below).




Each data file appears as a 2D representation of the run. At this stage you will be warned if any of the data files have been 'centroided' during the data acquisition and conversion process.

**Note:** the **Experiment Properties** are available from the File menu. These were selected when the experiment was created (see Appendix 1, page 54).



**Tip:** the **'Exclude areas from selected run'** facility allows you to examine and exclude areas (usually early and/or late in the LC dimension (Retention Time)) that appear excessively noisy due to capture of data during column regeneration (see Appendix 2, page 59). This is not required for this data set.



**Note:** use the **Remove Run** to remove run(s) from the current experiment.

Now start the Alignment process

## Stage 2: Automatic Alignment of your data

The process of alignment of your data can be started before the import of your data has been completed by clicking on **Start alignment process**.


**2 Start the alignment process**  
While your runs are importing, click the button below to:

- Select [alignment reference](#) candidates
- Determine the best of the candidates
- Align all runs to that reference run

Start alignment process

During this process the software will Align all your runs to a Reference run which can either be selected automatically by the software or manually selected by you.

In this tutorial example you have 6 runs which have already been imported so to start the process click on **Start alignment process**



**Progenesis LC-MS provides three methods** for choosing the alignment reference run, as seen below:

### 1. Assess all runs in the experiment for suitability

This method compares every run in your experiment to every other run for similarity.

The run with the greatest similarity to all other runs is chosen as the alignment reference.

If you have no prior knowledge about which of your runs would make a good reference, then this choice will normally produce a good alignment reference for you. This method can take a long time

### 2. Use the most suitable run from candidates that I select

This method asks you to choose a selection of reference candidates, and the automatic algorithm chooses the best reference from these runs.

When you have some prior knowledge of your runs suitability as references:

runs from pooled samples

runs for one of your experimental conditions will contain the largest set of common peptides.

### 3. Use this run


This method allows you to manually choose the reference run.

Manual selection gives you full control, but there are a couple of risks to note:

- If you choose a pending run which subsequently fails to load, alignment will not be performed.
- If you choose a run before it fully loads, and it turns out to have chromatography issues, alignment will be negatively affected (for this reason we recommend that you let your reference run fully load and assess its chromatography before loading further runs).


For this tutorial we will select the first option (See Appendix 1, page 54 for more details on using the other options).

You will now be asked if you want to Align your runs automatically or manually.



Select automatically and click finish.

The Alignment process starts with the automatic selection of C1 as the reference



Once the Reference run has been chosen the automatic alignment is then performed. As the whole process proceeds you get information on what stage has been performed and also the % of the process that has been completed.


When the Alignment completes you can either review the chromatography or go to the Review Alignment using the options on the Alignment Dialog.

Click **Review Alignment**.



### Stage 3: Licensing

This stage in the analysis workflow will **only** appear if you are using 'Unlicensed' data files to evaluate the software and have no dongle attached.



For details on how to use Licensing go to Appendix 3 (page 60)

**If you are using the tutorial archive, this page will not appear as the data files are licensed.**

### Stage 4: Review Alignment

At this stage Progenesis LC-MS Alignment opens displaying the alignment of the runs to the Reference run (C1).

**Table of Alignment Vectors and Scores**

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/> <input type="checkbox"/>	221	84.3%
A2	<input checked="" type="checkbox"/> <input type="checkbox"/>	240	82.2%
A3	<input checked="" type="checkbox"/> <input type="checkbox"/>	244	83.9%
C1	<input checked="" type="checkbox"/> <input type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/> <input type="checkbox"/>	362	97.8%
C3	<input checked="" type="checkbox"/> <input type="checkbox"/>	422	99.0%

### Layout of Alignment


To familiarize you with Progenesis LC-MS Alignment, this section describes the various graphical features used in the alignment of the LC-MS runs

To setup the display so that it looks similar to the one above:

- In the Run table click on Run A2 to make it current. You will now be looking at the alignment of A2 to C1 in the Unaligned view.

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/> <input type="checkbox"/>	221	84.3%
A2	<input checked="" type="checkbox"/> <input type="checkbox"/>	240	82.2%
A3	<input checked="" type="checkbox"/> <input type="checkbox"/>	244	83.9%
C1	<input checked="" type="checkbox"/> <input type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/> <input type="checkbox"/>	362	97.8%
C3	<input checked="" type="checkbox"/> <input type="checkbox"/>	422	99.0%

- Now adjust the size and position of the current focus. First select the size by clicking on the Focus grid size. Darken or lighten the runs using the contrast buttons. Then click on the Ion Intensity Map to 'locate' the current focus. The other 3 views will update to reflect the new focus.



**Vector Editing (Window A):** is the main alignment area and displays the area defined by the current **focus** rectangle shown in Window C. The current run is displayed in green and the chosen reference run is displayed in magenta. Here is where you can review in detail the vectors and also place the manual alignment vectors when required.


**Transition (Window B):** uses an **alpha blend** to animate between the current and reference runs. Before the runs are aligned, the features appear to move up and down. Once correctly aligned, they will appear to pulse. During the process of adding vectors, this view will change to show a zoomed view of the area being aligned to help accurate placement of manual vectors.

**Whole Run (Window C):** shows the **focus** for the other windows. When you click on the view the orange rectangle will move to the selected area. The focus can be moved systematically across the view using the cursor keys. The focus area size can be altered using the controls in the bottom left of the screen or by clicking and dragging out a new area with the mouse. This view also provides a visual quality metric for the Alignment of the runs (**note:** this can be switched off using the options in the View menu) which focuses your review of the alignment process.

**Total Ion Chromatograms (Window D):** shows the current **total ion** chromatogram (green) overlaid on the Reference chromatogram (magenta). As the features are aligned in the **Vector Editing** view the chromatograms become aligned. The retention time range displayed is the vertical dimension of the Focus Grid currently displayed in the **Whole Run** view (Window C).

This view assists in the verification of the feature alignment.

**Note:** the icon to the right of the 'Window' titles expands

Total Ion Chromatograms  the view .

## Reviewing quality of alignment vectors

After **Automatic alignment** the number of vectors and Quality Scores will be updated on the **Runs** panel and the vectors will appear (in blue) on the view.

If the alignment has worked well then in Windows A and C the grid lines (option under **View** menu) should show minimal distortion, Window B (Transition) will show features pulsing slightly but not moving up and down.

The screenshot shows the Progenesis LC-MS software interface. The 'Review Alignment' panel is active, displaying a table of runs and their alignment statistics. The 'Ion Intensity Map' is shown on the right, with a red box highlighting a specific area of interest. The interface includes a menu bar with options like 'File', 'Review Alignment', 'Filtering', etc., and a toolbar with buttons for 'Show Aligned', 'Show Unaligned', and 'Remove Vectors'.

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/>	221	84.3%
A2	<input checked="" type="checkbox"/>	240	82.2%
A3	<input checked="" type="checkbox"/>	244	83.9%
C1	<input checked="" type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/>	362	97.8%
C3	<input checked="" type="checkbox"/>	422	99.0%

Contrast: [Slider controls]

Ion maps: ■ Alignment target ■ Run being aligned

Focus grid size:  Make box square

Alignment quality: ■ Good ■ OK ■ Needs review



Section Complete

**Note:** you can use the icon to the right of the panel name to expand or contract each view

## Reviewing Quality of Alignment



At this point the quality metric, overlaid on the Ion Intensity Map as coloured squares, acts as a guide drawing your attention to areas of the alignment. These range from Good (Green) through OK (Yellow) to Needs review (Red). When reviewing individual squares set the grid size to 16, (and untick the Make box square option) using the '**Focus grid size**' control at the bottom left of the window. Three example squares are examined here.

For a 'green' square the majority of the data appears overlapped (black) indicating good alignment. When viewed in the Transition view the data appears to pulse.







For a 'yellow' square some of the data appears overlapped (black) indicating OK alignment. When viewed in the Transition view some of the data appears to pulse.



For a 'red' square little of the data appears overlapped (black) indicating questionable alignment. When viewed in the Transition view little data appears to pulse.



**Note:** the coloured metric **should be used as a guide**. In cases where there are a few 'isolated' red squares this can also be indicative of 'real' differences between the two runs being aligned and should be considered when examining the overall score and surrounding squares in the current alignment.

The weighted average of the individual squares gives the overall percentage score for each run alignment.

**Note:** a marked red area combined with a low score clearly indicates a 'miss alignment' and may require some manual intervention (see Appendix 4, page 61).

### Review Alignment

Align peptide ions to compensate for drifts in retention time by dragging them up or down in the Vector Editing window.

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/> <input type="checkbox"/>	1	14.6%
A2	<input checked="" type="checkbox"/> <input type="checkbox"/>	240	82.2%
A3	<input checked="" type="checkbox"/> <input type="checkbox"/>	244	83.9%
C1	<input checked="" type="checkbox"/> <input type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/> <input type="checkbox"/>	362	97.8%
C3	<input checked="" type="checkbox"/> <input type="checkbox"/>	422	99.0%

Contrast

Ion maps: ■ Alignment target ■ Run being aligned

Focus grid size:  Make box square

Show Aligned
Show Unaligned

#### Ion Intensity Map



## Stage 5A: Filtering

Now that you have reviewed your aligned Runs, you are ready to analyse them. Move to the **Filtering** stage, by either clicking on **Section Complete** (bottom right) or on Filtering on the workflow.



### Peak Picking Parameters

The Peak Picking Parameters dialog opens, showing all the runs in the current experiment and a tick against each run. This is the default setting, where the peak picking algorithm uses information from all of the runs to contribute to the pattern of feature outlines.




**Tip:** It may be appropriate **only** to pick peaks that are present in a limited number of your runs. In which case un-tick the runs that you do **NOT** want to contribute to the feature detection pattern. This may be important when one or more of the runs appear noisy due to non-optimal chromatography or sample handling.

**Note:** features outlines will be added to 'un-ticked' runs; although these runs will not contribute to the peak picking pattern.

**Tip:** depending on run quality, a suggested minimum number of ticked runs should include at least one replicate of each experimental condition.

The sensitivity of the peak picking can be controlled by adjusting settings under the **Peak picking limits** tab.

**Note:** for all 3 **Sensitivity** modes a Chromatographic peak width (Retention time window) for the peaks **can** be set by applying a minimum retention window or peak width in minutes.





**For the runs in this user guide we will use the default settings for the Automatic method and NOT apply a minimum peak width (as shown above).**

The third tab allows you to set the **maximum charge** of the peptide ions, which will be detected. The default setting is a charge state of 20. For this example leave this set as default.


**Note:** you can either, leave this set as default and remove the high charge state peptide ions at the Filtering stage or you can choose not to detect them in the first place by reducing the charge state threshold here.

Finally you can set **Retention time limits** for the detection. The default limits are displayed in the unticked boxes and detection will use these values whenever the box(s) are unticked. Enter values of 10 and 75 min and tick the boxes as shown below.



Click **Start peak picking** to start the detection process.

On completion of detection, the Filtering stage will open displaying the number of features detected, in this example there are 14624 features.



If required you can remove features based on position, charge state, number of isotopes or combinations of these feature properties.

As an example we will filter the features based on charge 'charge state'.

When **With charge** is selected the number of features present at each charge state is displayed, these can be selected accordingly. In this case we will retain features with a charge state of 2 to 7.

Area limits, charge state and number of isotopes can be combined to refine the feature selection.

**Tip:** when filtering on one property of the feature i.e. charge state, make sure you have 'collapsed' the other filters (see right)

**Filter Features**

You can filter features that you do not wish to include in your analysis by using the criteria below.

Select all features matching the following filters:

- Inside area
- With charge
  - Charge 1 (911 features)
  - Charge 2 (5904 features)
  - Charge 3 (5357 features)
  - Charge 4 (1729 features)
  - Charge 5 (483 features)
  - Charge 6 (73 features)
  - Charge 7 (46 features)
  - Charge 8 (29 features)
  - Charge 9 (22 features)

Therefore all features with a charge state of 1 or 8 and above will appear blue on the main view.




To remove these features press **Delete 1032 Non Matching Features**.

You can use the **Undo** button to bring back deleted features, however, when you move to the next section you will lose the capacity to undo the filter. Before moving on from filtering you can review the normalisation of the data.


**Tip:** When you have reached the filtering stage, it is good practice to close the experiment and save an archive. This can be used to restore the unfiltered state if the filtering you have performed is too 'stringent'.

## Stage 5B: Reviewing Normalisation

Normalisation review is accessed from the button at the bottom right corner of the filtering page.



If you have filtered out a number of features from the original detection pattern then the normalisation will update.




The **Review Normalisation** page will open displaying plots for the normalisation of all the features on each run.

This page in the workflow **does not** allow you to alter the Normalisation of your data but provides you with individual views for each run showing the data points used in the calculation of the normalisation factor for the run.

Alternatively, if you do not believe normalisation is necessary, you can opt to 'not use any normalisation' for the rest of the analysis (Normalisation Method tab).




Normalisation factors are reported in the table to the left of the plots.



**Calculation of Normalisation Factor:**


Progenesis LC-MS will automatically select one of the runs that is 'least different' from all the other runs in the data set to be the 'Normalising reference'. The run used, is shown above the table of Normalisation factors.

For each sample run, each blue dot shows the log of the abundance ratio for a different feature (normalisation target abundance/run abundance).




The details for individual features can be viewed as you hold the cursor over the dots on the plot.

On the graph the features are shown ordered by ascending mean abundance. The normalisation factor is then calculated by finding the mean of the log abundance ratios of the features that fall within the 'robust estimated limits' (dotted red lines). Features outside these limits are considered to be outliers and therefore will not affect the normalisation.




Finally, if you do **not** wish to work with normalised data then Select '**Don't use any normalisation**'.



**Note:** once you have identified a set of proteins, you can then apply the **Normalise to a set of housekeeping proteins** by using this option to locate and select the features.

For this experiment, you should leave the **Normalise to all features** option selected.

Now return to filtering by clicking on the button on the bottom left of the screen



For this example, we **DO NOT** do any additional Filtering so click on **Section complete**.


**Note:** if you do any extra filtering then **Normalisation recalculates** as you move to the next stage in the Workflow.

## Stage 6: Experiment Design Setup for Analysed Runs

At this stage in the workflow you can setup one or more experimental designs for your data.

There are two basic types of experimental designs:

**Between-subject design:** here samples from any given subject appear in only one condition. (i.e. control versus various drug treatments). The ANOVA calculation assumes that the conditions are independent and applies the statistical test that assumes the means of the conditions are equal.



LC-MS Tutorial - Progenesis LC-MS

File Import Data Review Alignment Filtering Experiment Design Setup Review Peak Picking Peptide Statistics Identify Peptides Refine Identifications Resolve Conflicts Review Proteins Protein Statistics Report nonlinear DYNAMICS Help

New

Which experiment design type do you want to use for this experiment?

**Between-subject Design** Create

Do samples from a given subject appear in only one condition? Then use the between-subject design.

To set up this design, you simply group the runs according to the condition (factor level) of the samples. The ANOVA calculation assumes that the conditions are independent and therefore gives a statistical test of whether the means of the conditions are all equal.

	A	Delete
	A1	Remove
	A2	Remove
	A3	Remove
	C	Delete
	C1	Remove
	C2	Remove
	C3	Remove
	<a href="#">Add condition...</a>	

**Within-subject Design**

Have you taken samples from a given subject under different conditions? Then use the within-subject design.

**Note:** you must have a sample from every subject for every condition to use a within-subject design.

For example, you would choose this type of design for a time series experiment where every subject has been sampled at each time point.

To set up this design, you tell the software not only which condition (factor level) each run belongs to but also which subject it came from. The software will then perform a repeated measures ANOVA.

A standard ANOVA is not appropriate because the data violates the ANOVA assumption of independence. With a repeated measures ANOVA individual differences can be eliminated or reduced as a source of between condition differences (which helps to create a more powerful test).

The within-subject design can be thought of as an extension of the paired-samples t-test to include comparison between more than two repeated measures.


	Before	During	After
Patient X	X1	X2	X3
Patient Y	Y1	Y2	Y3
Patient Z	Z1	Z2	Z3

**Within-subject design:** here samples have been taken from a given subject under different conditions (i.e. the same subject has been sampled over a period of time or after one or more treatments). Here a standard ANOVA is not appropriate as the data violates the ANOVA assumption of independence. Therefore by using a *repeated measures* ANOVA, individual differences can be eliminated or reduced as a source of between condition differences. This within-subject design can be thought of as an extension of the paired samples t-test, including comparison between more than two repeated measures.


Additional information on how to apply the Within-subject Design is in **Appendix 5** page 66

This experiment contains 2 conditions: A and C and uses the **Between-subject design** to group the analysed runs to reflect the Biological conditions in the original study.

To create a new **Between-subject Design** move the cursor over this option and click to open the dialog.




Give the new experimental design a name and then click **Create design**.



## To create a new condition


1. Select the runs for the condition by clicking on the required icon in the **Runs** panel, as shown.
2. Press the 'black triangle' next to the **Add Selected Runs to Condition** button on the main toolbar.
3. Select **Add to new condition...** from the drop down menu.
4. A new condition will appear in the **Conditions** panel on the left.
5. Rename the condition (e.g. C) by over typing the default name.
6. Repeat steps 1 to 5 until all the required runs are grouped into conditions.

An alternative way to handling the grouping of this set and other larger (and more complex) experimental designs is to make use of **sample tracking information** that has been stored in a spreadsheet at the time of sample collection and/or preparation.




For this example there is a **Tutorial Groups.csv** file available in the Experiment Archive you restored at the beginning of this tutorial exercise.

To use this approach select the **Import design from file** option from the New Experiment Design dialog. Then locate the Tutorial Groups file and select what to **Group by**, for example: **Conditions**.



When **Create design** is pressed the new tab refreshes to allow you to adjust the conditions.

Use **Delete** on the Conditions panel to remove conditions that are not required in this particular design.



**Note:** On deleting a condition the runs will reappear in the Runs window.

**Note:** both designs are available as separate tabs.

To move to the next stage in the workflow, **Review Peak Picking**, click **Section Complete**.

## Stage 7: Review Peak Picking and editing of results

The purpose of this stage in the Workflow is to review the list of features using the visual tools provided and edit features if required.

The review stage has 4 display modes: 1D, 2D, 3D and Feature Details controlled by the tabs on the bottom left of the display and the expander bar to the right of the table. Each display has multiple views to allow comparative exploration of the detected features on the aligned LC-MS runs.

### Exploring analysed data using the Data displays

**Window A:** shows the list of features ranked by the p value for the one way **Anova** using the current grouping.

**Note:** a value of 'Infinity' in the **Fold** column indicates 'Presence/Absence'

The screenshot displays the Progenesis LC-MS software interface. The main window is titled "Review Features" and shows a table of features ranked by p-value. The table has columns for #, Anova (p), Fold, Tag, and Notes. The feature at row 1391 is highlighted with a red box. To the right of the table are two plots: a Mass spectrum (25.179 min) and a Chromatogram (m/z = 1099.5684). Below the table are buttons for "Delete 1 selected feature" and "1D Display", "2D Montage", "3D Montage". At the bottom right, there is a "Section Complete" button.

#	Anova (p)	Fold	Tag	Notes
6314	6.47E-11	Infinity		Adc
10970	6.75E-11	Infinity		Adc
11767	7.27E-11	Infinity		Adc
1455	9.01E-11	3.5E+08		Adc
8761	1E-10	Infinity		Adc
8045	1.01E-10	Infinity		Adc
1391	1.01E-10	Infinity		Adc
1142	1.21E-10	Infinity		Adc
14023	1.27E-10	Infinity		Adc
11985	1.36E-10	Infinity		Adc
7057	1.48E-10	Infinity		Adc
8028	1.56E-10	Infinity		Adc


**Note:** by default all the features are included in the selection for the next section of the analysis.

To highlight a group of features drag out a selection on the table.




**The 1D Display**

**Window B:** displays the Mass spectrum for the current feature on the selected Run (in window D).




**Window C:** displays the Chromatogram for the current feature on the selected Run (in window D).




**Window D:** displays the details of the currently selected run. By default the selected run is an Aggregate of all the aligned runs.

Details of individual runs can be viewed by using the 'Run' link and selecting the run you wish to view.




The feature editing tools are located in this window (see page 26 for functional explanation).

Clicking on the Expression Profile tab in Window D shows the comparative behaviour of the feature across the various biological groups based on group average normalised volume. The error bars show +/- 3 standard errors.




**Window E:** shows where the current feature is located on the LC-MS run by means of the 'Green' rectangle.

To change the current location, click on the image of the run (**note:** the retention time and m/z values update as you move the cursor around this view).



**Note:** doing this updates the focus of all the other windows.


You can also drag out a larger area on this view that will refocus the other windows.



### The 2D Display

Windows **A**, **D** and **E**: perform the same functions across all 4 display modes.


In the 2D Montage mode, Window B displays a montage of the current feature across all the aligned LC-MS runs.



The appearance of the Montage (window B) is controlled by the panel on the bottom left of the display.




Using the the various views in the 2D display one can examine the feature detection in detail to validate the correct detection of even fully overlapping features as shown above.



### The 3D Display


Window B changes into a 3D view by selecting the 3D Montage tab on the bottom left of the display.




The number of 3D views displayed in the montage is controlled using the [Select runs](#) link on the 3D Montage tab. The views can be set to **Rotate** automatically or you can rotate them manually by clicking and dragging them with the mouse.

## Editing of features in the View Results stage


As an example of using the editing tools which are located on the left of the LC-MS Run view, we will remove and add back the 'monoisotopic peak' for the detected feature selected below. A feature can be selected from the 'Features' list or located using the various views.




1. Locate the feature at approx 1004.77 m/z and 30.453 min using the **Find** tool.



2. Select the **Edit** tool and click on the feature to reveal the 'edit handles'




3. Click on the 'minus' handle over the monoisotopic peak to remove it.




4. Click outside the boundary of the feature to update the view.

- To add a peak to an existing feature, ensure that **Edit** is selected then click inside the feature to reveal the handles.



- Click on the 'plus' handle on the peak to add it.
- Then click outside the feature to update the view.
- Note:** If you are not satisfied with the editing use the **Undo** button and retry.
- Finally note: that a tag is automatically added to the edited feature in the table and the features id number is changed to the next available one at the end of the list.




The other tools: split, merge, add and delete behave in a similar fashion and their use can be combined to achieve the desired results.

### Selecting and tagging features for Peptide Statistics

There are a number of ways to 'refine' your 'Ranked List' of analysed features before examining them with the Statistical tools in **Peptide Statistics**. These make use of simple 'Selection' and 'Tagging' tools that can be applied to the various groupings created in Stage 6 (page 19). An example is described below.

First expand the 'Features' table to show all the details by clicking on the 'Expander bar' to the right of the Review Features table.



Then order on **Abundance** and select all features with an Abundance of 1E+05 and greater, (the exact number is not important).

The screenshot shows the 'Review Features' window in Progenesis LC-MS. The 'Abundance' column is highlighted in the table header. The table contains the following data:


#	Anova (p)	Fold	Tag	Notes	Highest	Lowest	m/z	z	Mass	RT (mins)	RT window	Abundance	Intensity	Max CV
2	0.159	1.09		Add note...	A	C	805.441	3	2413.301	54.875	7.65	1.25E+08	1.08E+08	8.15
18	0.463	1.05		Add note...	C	A	1207.6552	2	2413.296	54.899	4.9	9.09E+07	5.04E+07	8.52
24	3.18E-06	5.72E+04		Add note...	C	A	1100.5867	3	3298.738	44.928	2.43	8.39E+07	7.87E+07	53.4
78	3.41E-06	1.15E+03		Add note...	A	C	1176.227	3	3525.659	48.238	9.44	6.91E+07	2.22E+07	23.1
14	4.37E-07	856		Add note...	A	C	656.8612	2	1311.708	44.137	4.21	6.19E+07	1.17E+08	15
20	2.22E-07	786		Add note...	C	A	988.9849	2	1975.955	50.605	4.65	5.14E+07	9.19E+07	14.4
7	8.31E-08	130		Add note...	C	A	663.8693	2	1325.724	46.597	4.14	5E+07	1.69E+08	6.68
23	2.08E-06	5.15E+03		Add note...	C	A	900.9713	2	1799.928	39.272	3.02	4.57E+07	8.16E+07	37.5
56	5.85E-06	3.79E+03		Add note...	A	C	1061.007	2	2119.999	53.288	6.75	4.3E+07	2.52E+07	33.6
41	5.57E-06	1.33E+03		Add note...	A	C	997.4477	2	1992.881	31.814	2.57	4.24E+07	3.81E+07	34.3
94	0.0589	2.84		Add note...	C	A	976.8129	3	2927.417	54.357	13.3	3.72E+07	1.58E+07	55.3
47	7.38E-06	3.58E+03		Add note...	C	A	1032.4669	3	3094.379	32.959	3.63	3.53E+07	3.95E+07	32.8

With the features still highlighted right click on them and select 'New Tag'.

The screenshot shows the 'Review Features' window with a context menu open over a feature. The 'New tag...' option is highlighted. The table shows the following data:


#	Anova (p)	Fold	Tag	Notes	Highest	Lowest	m/z	z	Mass	RT (mins)	RT window	Abundance	Intensity	Max CV
2929	0.00372	6.81		Add note...	C	A	630.866	5	3149.294	33.122	1.23	1.01E+05	4.42E+05	40.8
1833	0.0187	3.03		Add note...	A	C	681.8701	2	1361.726	52.984	0.666	1.01E+05	4.15E+05	42.8
4035	0.0203	1.74		Add note...	C	A	1196.9362	3	3587.787	39.3	0.921	1E+05	2.57E+05	20.8
5464	2.85E-08	Infinity		Add note...	C	A	1585.5178	3	4753.532	46.597	0.359	1E+05	1.7E+05	16.9
1686	0.0798	2.62		Add note...			706.6327	4	2822.502	63.027	1.26	1E+05	5.39E+05	62.6
4771	0.0177	540		Add note...			1157.6035	3	3469.789	46.561	0.808	1E+05	2.24E+05	126
2941	0.0748	1.49		Add note...			1081.0327	2	2160.051	36.226	0.754	1E+05	1.96E+05	25.4
3391	0.000569	3.03		Add note...			703.845	2	1405.675	45.667	1.94	9.99E+04	2.53E+05	17.7
5681	0.00365	20.8		Add note...			1019.4452	3	3055.314	53.762	0.925	9.99E+04	1.94E+05	69.3
8312	0.000149	6.8		Add note...			643.3413	2	1284.668	18.939	1.32	9.99E+04	2.37E+06	16.9
2381	0.00089	4.49		Add note...	A	C	817.9225	2	1633.83	39.467	0.743	9.98E+04	2.87E+05	22.4
2379	0.0448	150		Add note...	C	A	814.8726	4	3255.461	37.56	0.71	9.98E+04	4.7E+05	118

Give the Tag a name. i.e. 'Most abundant'.




On clicking **OK** the Tag is added to the features highlighted in the table (signified by a coloured square).

Now right click on any feature in the table and select **Quick Tags** this will offer you a number of standard tag options. Select **Anova p-value....** Then set the threshold as required and adjust default name as required and click **Create Tag**.




Once this tag appears against features in the table right click on the table again and create another Quick Tag, this time for features with a **Max fold change ≥ 2**






The table now displays features with multiple tags.

The tags can be used to quickly focus the table on those features that display similar properties.

For example: to focus the table on displaying those features that have an **Anova p-value  $\leq 0.05$**  click on **Create** on the filter panel above the table.




Drag the tag on to the panel **Show features that have all of these tags** and press OK.



Now order the current features in the table by the **Highest mean** so that all the features showing the highest mean for **condition C** are at the top of the list.

Then highlight all the features with the highest mean for **condition C** and create a new Tag for them.


#	Anova (p)	Fold	Tag	Notes	Highest	Lowest
1440	2.94E-05	74.5		Add note...	C	A
951	0.00113	138		Add note...	C	A
623	0.000379	12.6		Add note...	C	A
3206	0.0224	11.3		Add note...	C	A
643	1.73E-05	112		Add note...	C	A
692	0.000628	411		Add note...	C	A
853	7.89E-06	5.18E+05		Add note...	C	A
11868	0.0289	49.9		Add note...	A	C
5801	0.000122	21.3		Add note...	A	C
13055	9.77E-08	Infinity		Add note...	A	C
1842	0.00615	2.99		Add note...	A	C
5030	0.000154	32.2		Add note...	A	C




#	Anova (p)	Fold	Tag	Notes	Highest	Lowest
1440	2.94E-05	74.5		Add note...	C	A
951	0.00113	138		Add note...	C	A
623	0.000379	12.6		Add note...	C	A
3206	0.0224	11.3	Most abundant	Add note...	C	A
643	1.73E-05	112	Anova p-value <= 0.05	Add note...	C	A
692	0.000628	411	Max fold change >= 2	Add note...	C	A
853	7.89E-06	5.18E+05		Add note...	C	A
11868	0.0289	49.9		Add note...	A	C
5801	0.000122	21.3		Add note...	A	C
13055	9.77E-08	Infinity		Add note...	A	C
1842	0.00615	2.99		Add note...	A	C
5030	0.000154	32.2		Add note...	A	C

Create a tag for them called **Significantly up in C**, tagging 3959 features.

Finally view the tags you have just created by clicking on Edit in the Tag filter panel, above the table.



Make sure that only the tag for the **Most abundant** features is shown and press **OK**.



To move to the next stage in the workflow, Peptide Statistics, click **Section Complete**.

## Stage 8: Peptide Statistics on selected features


The user guide now describes the functionality of the Multivariate Statistics.



Peptide Statistics opens calculating the Principal Components Analysis (PCA) for the active 'tag' in this case the **Most abundant** features.

As an example we will start by examining the behaviour of the **Most abundant** features from the previous stage, **Review Peak Picking**.

The statistical analysis of the selected data is presented to you in the form of interactive graphical representation of answers to questions asked of the analysed data.



**Note:** the LC-MS runs (samples) are displayed as solid coloured circles on the plot. To identify the runs, a tooltip is displayed when the cursor is held over each circle.

### Principal Component Analysis (PCA)

In **Peptide Stats** the first statistically based question asked of the data takes the form of a Quality Control assessment:


*Are there any outliers in my data? And does my data cluster according to my experimental conditions?*

It answers this question by:



'Using Principal Components Analysis (PCA) to produce a simplified graphical representation of your multidimensional data'.

PCA can be used to determine whether there are any outliers in the data and also look at how well the samples group. The groupings that can be observed on the 2D PCA plot can be compared to your experimental conditions and conclusions can be drawn regarding possible outliers in your data. Selecting features in the table will highlight the features on the 'Biplot' and their expression profiles will appear in the lower panel.



**Note:** the Table in the Stats view contains additional columns:

**q value:** tells us the expected proportion of false positives if that feature's p-value is chosen as the significance threshold.

**Power:** can be defined as the probability of finding a real difference if it exists. 80% or 0.8 is considered an acceptable value for power. The Power Analysis is performed independently for each feature, using the expression variance, sample size and difference between the means.

Also, for a given power of 80% we can determine how many samples are required to ensure we find a difference if it actually exists.

**Note:** Power analysis is discussed in Appendix 6 (page 68)




## Correlation Analysis

Use the tags, created in Review Peak Picking, to filter the features in the table.


We are going to explore the Correlation Analysis for all the features that were tagged at the view results stage for having an **Anova p-value ≤ 0.05**.

On pressing OK the PCA will recalculate using these 8618 features, you can (to save time) stop this calculation by pressing **Cancel calculation** and then set up Correlation Analysis for the 7965 features.



To set up the **Correlation Analysis** using this filtered data set click on **Ask another question** (above the table)

A selection of 3 tools will appear in the form of questions.



Select the second option to explore ‘feature correlation based on similarity of expression profiles’




This time the statistically based question(s) being asked is:

*‘Group my (selected) features according to how similar their expression profiles are’*

The question is answered by:


*‘Using Correlation analysis to evaluate the relationships between the (selected) features’ expression profiles’.*

The answer is displayed graphically in the form of an interactive dendrogram where the vertical distance, between each feature can be taken as indicative of how similar the expression profiles of each cluster of features are to each other.




**Correlation Analysis** enables the grouping of features together according to how similar their expression profiles are.

To highlight all the features demonstrating **Increased expression in the C** group click on a 'node' for a branch of the Dendrogram (as shown above). As before right click on the highlighted features in the table and create a Tag for these features (Up regulated in C).




Also create a tag for those features showing **Increased expression in A** by first clicking on the other 'main' node then right click on the highlighted features in the table and creating the New tag (Up regulated in A).



**Comment:** When you review the tags using **Edit** you can see that the Magenta and Yellow tags have been assigned to the same number of features. This shows how tabulated information about features can be used alongside interactive graphical plots of multivariate statistical analysis to explore your data.

**Note:** two groups is a special case, for more groups this will not be the case additional tagging will be required.




**Tip:** when reviewing the tags (see above) if you are not applying a new filter then use the **Cancel** button to return to the main view, this prevents unnecessary recalculation of your data.

To move to the next stage in the workflow, **Identify Peptides**, click **Section Complete**.

## Stage 9: Identify peptides

Progenesis LC-MS does not perform peptide identifications itself. Instead it supports identifications by allowing you to export MS/MS peak lists in formats which can be used to perform peptide searches by various search engines. The resulting identifications can then be imported back into Progenesis LC-MS, using a number of different file types, and matched to your detected features.



Determining protein identification is dependent on the availability of MS/MS data for the LC-MS runs. This data may be available but limited if the LC-MS was performed in a data dependant MS/MS detection mode due to under sampling. Under these conditions MS/MS data acquisition is dependent on thresholds and parameters set prior to the acquisition of the LC-MS run.

For this example we are using LC-MS runs containing MS/MS data where the data was acquired in a data dependant mode (DDA).

The Peptide Search page shows the number MS/MS that have been matched to each feature in the Feature list (see above). MS/MS scans are matched to a feature if their precursor *m/z* and aligned retention time fall within the area of one of the isotopes of the feature. The MS/MS scans which are matched to the displayed features are shown in the MS/MS spectra list on the right.

The first step is to decide which MS/MS scans you wish to export to be identified. By default this is all the available spectra for the features displayed in the Features list (in this case all the features that have an **Anova p-value ≤ 0.05**. This number is visible on the Export button.

The set can be targeted using the tags and also refined with respect to quantity and quality of the spectra being sent to the search engine.

Filter the table to show only the features tagged **Significant p<0.05** as shown.

**Note:** by default the table is ordered on the number of MS/MS spectra available for each feature.

The total number of spectra included in this set is **20769** as shown on the Export button.

Before exporting the spectra, the set can be further refined.

**Note:** many of the abundant features have a large number of spectra associated with them.

**Performing the search**

Select the search program you're using:  
 Mascot

[Help](#)

**Export 20769 ms/ms spectra**

MSMS Preprocessing

Limit fragment ion count 40

Deisotoping and charge deconvolution

Import search results

Clear all identifications

**Identify Peptides**

Run ms/ms ion searches by exporting peak list files to a protein identification program. Search results files may then be imported back in to link protein IDs to detected peaks.

Features:

Tag filter applied: features may be hidden

#	MS/MS	Proteins	Score	Tags
56	50	0		
34	58	0		
7	54	0		
127	54	0		
78	52	0		
91	50	0		
64	46	0		
3820	45	0		
59	45	0		
135	43	0		
20	41	0		
21	41	0		
30	40	0		
32	39	0		
425	39	0		

**MS/MS Spectra**

Batch inclusion options for creating export list of ms/ms spectra

Rank: greater than 5

Feature ID: less than

Charge: less than

Scan number: less than

Exported: equal to

Isotope: less than

ID score: less than

Feature intensity: less than

Precursor intensity: less than

Precursor intensity (%): less than

Run name: contains

Peptide sequence: contains

Protein accession: contains

Protein description: contains

Include in export Exclude from export Clear all filters

Export	Rank	#	Run	Scan number	Exported	Feature intensity	Precursor intensity (%)	Charge	Precursor m/z	Isotope	Id
<input checked="" type="checkbox"/>	4	132	A3	4540	Yes	9.0e+006	3.8e+006	41.7	2	980.9709	1
<input checked="" type="checkbox"/>	5	132	A2	4602	Yes	1.4e+007	5.5e+006	39.8	2	980.9717	1
<input type="checkbox"/>	6	132	A1	4461	No	1.1e+007	3.7e+006	34.8	2	980.9718	1
<input type="checkbox"/>	7	132	A3	4580	No	9.0e+006	1.5e+006	16.3	2	980.9711	1

Feature number 129, m/z 952.7989, retention time 53.245 min, charge +3

Run: C1 Scan number: 6584

Retention time (min) vs m/z plot showing peaks at 953 and 954.

Intensity vs m/z plot showing a major peak at m/z 952.7989.

Section Complete

To control the number of spectra for each feature, expand the **Batch inclusion options**.

**For example:** We will make use of the 'Rank' value to reduce the number of Spectra being used for each feature to a maximum of 5.

The 'Rank' of each MS/MS spectra is determined by comparing its **% value against all other spectra matched to the same feature**.

Export	Rank	#	Run	Scan number	Exported	Feature intensity	Precursor intensity	(%)	Charge	Precursor m/z	Isotope	Id score
<input type="checkbox"/>	23	9	A1	4868	No	1.2e+008	1.1e+006	1.0	2	656.8616	1	
<input type="checkbox"/>	24	9	A2	5311	No	1.4e+008	1.4e+006	1.0	2	656.8618	1	
<input type="checkbox"/>	25	9	A3	5042	No	9.3e+007	8.8e+005	1.0	2	656.8614	1	
<input type="checkbox"/>	26	9	A1	5379	No	1.2e+008	1.0e+006	0.9	2	656.8610	1	
<input type="checkbox"/>	27	9	A3	5374	No	9.3e+007	7.7e+005	0.8	2	656.8615	1	

The rank of each MS/MS spectrum found by comparing its '% values against all other spectra matched to the same feature.

**Note:** the % value for each spectra is the **Precursor intensity as a percentage of the Feature intensity**

Set the Rank filter to 'greater than' 5 and click **Exclude from export** this reduces the number to spectra to export to **11544**

Limiting the 'fragment ion count' (FIC) for the spectra being exported can improve the quality of the spectral data being used in the search by removing noisy peaks.

For example for the current spectra, reduce the FIC from 1000 to 40.

The image shows two screenshots of the Progenesis LC-MS software interface. Both screenshots show the 'Performing the search' panel on the left and a spectral plot on the right. The top screenshot shows the 'Limit fragment ion count' set to 1000, and the bottom screenshot shows it set to 40. The spectral plots show a significant reduction in the number of peaks when the FIC is limited to 40. A red arrow points from the top spectrum to the bottom one, indicating the effect of the change.

**Note:** the effect this has on the number of peaks in the spectra. This 'limitation' is applied to all the spectra being exported; hence the export file size will be reduced.

For this example we will **not** limit the fragment count, so leave it un-ticked (the default setting).

## Performing an MS/MS Ion Search



Having chosen 11544 spectra to export, as described above:

1. Select appropriate search engine i.e. Mascot
2. Click 'Export current query set' to save search as file
3. Perform search on appropriate search engine and save results file
4. Click 'Import search results', locate results file and open

Please refer to Appendix 7 (pages 69) for details of the 'Search Engine' parameters used in this example

**Note:** the blue link provides you with details on the appropriate formats for exporting search results and access to additional formats

**Note:** an example Search Results file, from a MS/MS Ion search, is available in the folder you restored the Archive to (Protein Search Results.xml). Select the 'Mascot' method and import this file to see results like those below.



On importing the Search results the Features table updates to reflect the identified proteins and the relevant score for each searched feature.

**Identify Peptides**  
Run [ms/ms ion searches](#) by exporting peak list files to a protein identification program. Search results files may then be imported back in to link protein IDs to detected peaks.

Features:

#	MS/MS	Proteins	Score	Tags
56	59	1 gjj2170250...	102	
34	58	2 gjj2549763...	37.7	
7	54	6 gjj1121811...	32.8	
127	54	0		
78	52	3 gjj2551019...	111	
91	50	2 gjj2549763...	72.5	
64	46	0		
3820	45	0		
59	45	0		
135	43	4 gjj2170250...	32.3	
20	41	5 gjj2067250...	34.2	
21	41	1 gjj3661641...	44.5	
30	40	3 gjj2170250...	30.9	
32	39	4 gjj2170250...	57.5	
425	39	0		

**MS/MS Spectra**  
Batch inclusion options for creating export list of ms/ms spectra

Rank: greater than 5  
Feature ID: less than  
Charge: less than  
Scan number: less than  
Exported: equal to  
Isotope: less than

Feature intensity: less than  
Precursor intensity: less than  
Precursor intensity (%): less than  
Run name: contains  
Peptide sequence: contains  
Protein accession: contains  
Protein description: contains

Include in export Exclude from export Clear all filters

**Import search results**  
6978 search hits have been imported and assigned to features

Precursor intensity (%)	Charge	Precursor m/z	Isotope	Id
6.2e+007	83.5	2	595.3201	1
9.4e+007	81.7	2	595.3196	1
7.3e+007	66.7	2	595.3195	1

**Feature number 10, m/z 595.319, retention time 36.921 min, charge +2**

Run:A3 Scan number:3970

In order to review, and refine the quality of the **Peptide Search** results click on the next stage in the workflow, **Refine Identifications**.



## Stage 10: Refine Identifications

In this example the organism under study is *Clostridium difficile*

**Note:** before removing any identifications, make sure there are **no** tag filters applied at the Identify peptides stage

As an **example** 'Acceptance Criteria' on which to base the sequential filtering of the Peptide results, the following thresholds will be applied:

- Remove identifications with a Score less than 40
- Remove identifications where less than 2 hits were returned
- Remove all identifications where the Protein Description **Contains** 'hypothetical'
- Remove all identifications where the Protein Description **Doesn't contain** '*Clostridium difficile*'

The screenshot displays the 'Peptide Search Results' window in Progenesis LC-MS. The 'Batch deletion options' panel is active, with the 'Score' filter set to 'less than 40'. Other filters include 'Hits' (less than), 'Mass error (ppm)' (less than), 'M/Z' (less than), 'Retention Time' (less than), 'Sequence Length' (less than), 'Charge' (less than), 'Sequence' (contains), 'Accession' (contains), 'Description' (doesn't contain), and 'Modifications' (contains). The 'Delete matching search results' button is highlighted. The main table shows search results with columns for #, Score, Hits, m/z, RT(mins), Charge, Mass, Mass err, Sequence, Accession, and Modifications. A red box highlights the bottom of the table showing '1846 search results, 486 matching batch delete options.' The 'Section Complete' button is visible at the bottom right.


On the Batch detection options panel, set the Score to less than 40, then **Delete matching search results**.

**Note:** the search results matching the filter criteria turn pink and the total is displayed at the bottom of the table (486 matching out of 1846)

**Note:** a dialog warns you of what you are about to delete

Now **Clear all filters** and then apply the next filter (Hits: less than 2) followed by the remaining two filters (page 39)

Having applied all 4 filters the **Peptide Search Results** should be reduced to **1129**.



To validate the Peptide search results at the protein level select the next stage in the workflow by clicking on **Resolve Conflicts**.

## Stage 11: Resolve Conflicts

This stage allows you examine the behaviour of the identified peptides and resolve any conflicts for the various peptide assignments at the protein level.

The **Resolve Conflicts** stage provides a number of interrelated graphical and tabular views to assist you in the validation of the peptides that have been assigned to proteins and also to review the relevance of the data returned from the search.

Open **Resolve Conflicts** and order the data in the Proteins table (A) on the basis of **Conflicts**.

**Note:** the look of the tables (with regards to ordering) in the following section may vary slightly.

The screenshot shows the Progenesis LC-MS software interface. The 'Resolve Conflicts' stage is active. The 'Proteins' table (A) is sorted by 'Conflict'.

Accession	Peptides	Conflict	Score	Tags
gj 254976387 (+6)	33 (18)	27	3.75E+03	
gj 255101963 (+4)	30 (15)	27	3.32E+03	
gj 254976385 (+1)	11 (2)	17	1.06E+03	
gj 255656776 (+1)	9 (1)	17	786	
gj 126700407	9 (1)	16	886	
gj 209571234 (+3)	22 (12)	11	2.25E+03	
gj 384359782 (+8)	21 (11)	11	1.88E+03	

The 'Peptides of gj|254976387' table (B) shows the following data:

#	Score	Hits	Mass	Mass error (p...)	RT (mins)	Charge	Tags	Abundance	Conflict	Peptide S
607	85.7	5	2926.431	0.0963	52.1	4		6.56E+05	1	TYN...
91	72.5	5	2926.432	0.127	52	3		1.09E+07	1	TYN...
5672	102	4	1800.888	0.1	19	2		2.88E+04	1	VEI...
12	73.1	5	1144.634	0.1	41.7	2		2.28E+07	1	GLL...
16	73.1	5	1144.634	0.371	40.7	2		2.69E+07	1	GLL...
1097	67.8	5	1144.633	-0.576	38	2		2.69E+05	1	GLL...
566	125	5	1832.936	-0.565	24.3	2		9.91E+05	1	ILN...

The 'Protein Resolution' graph (C) shows 'Standardised Expression Profiles' for two conditions, A and C. The y-axis is 'Standardised Normalised Abundance' ranging from -2.0 to 1.5. The x-axis represents different peptides. The graph shows a significant increase in abundance for several peptides in condition C compared to condition A.


The 'Retention time (min)' vs 'm/z' plot (D) shows a zoomed-in view of the retention time range from 18.537 to 19.43 minutes and m/z from 901 to 903. A red box highlights a specific peak at approximately 19.0 minutes and m/z 902.

The 'Retention time (min)' vs 'm/z' plot (E) shows a scatter plot of retention time (min) versus m/z. The y-axis ranges from 0.01 to 86.966 minutes, and the x-axis ranges from 500 to 2000 m/z. A red box highlights a cluster of points around 41.695 minutes and m/z 900.

Depending on the ordering, make 'cell surface protein precursor etc' the current protein by clicking on it in Window A (a circular orange symbol indicates current protein).




This protein has 33 peptides assigned (window B) which have a total of 27 conflicts. To view the conflicting assignments click on the **Protein Resolution** tab (window C).



In this case the conflicting peptide assignments are with **'The same protein'** (from a different strain) which also contains 27 conflicts. A simple resolution to these conflicts is to right click on the conflicting protein and turn off all its peptides (based on lower number of peptides and score).

**Note:** as you un-assign the peptides the number of conflicts update 'on the fly' in all the windows.



In this case the conflicting peptides are unassigned from the 'precursor' protein.

In some cases you can resolve the conflicts between 2 proteins on the basis of consistent peptide expression. In the example below the proteins share 11 conflicts. For the protein showing 2 clear patterns of expression you can un-tick all the peptides with conflicts in the corresponding peptides table

LC-MS Tutorial - Progenesis LC-MS

Experiment design: AC

Proteins: No filter applied

Accession	Peptides	Conflict	Score	Tags
gi 384359782 (+8)	21 (11)	11	1.88E+03	
gi 209571234 (+3)	22 (12)	11	2.25E+03	
gi 126699078 (+2)	3 (1)	3	307	
gi 126699128	3 (1)	3	328	
gi 254976383 (+1)	5	0	428	
gi 255656774 (+2)	1	0	136	

Peptides of gi|209571234

#	Score	Hits	Mass	Mass error (p...)	RT (mins)	Charge	Tags	Abundance	Conflict	Peptide
2316	133	5	2446.223	0.0965	62.3	2		1.56E+05	1	NA
204	98.9	5	1502.794	-0.167	28.9	2		8.17E+05	1	GE
225	41.7	5	1732.899	-0.452	29.1	3		6.66E+05	1	A.T.
461	77.3	5	1732.899	-0.502	29.1	2		5.09E+05	1	A.T.
288	52.3	5	1051.555	-0.192	23.4	2		4.82E+05	1	Y.Q.
334	127	5	1692.929	-0.187	45.4	2		8.22E+05	1	V.Y.

Protein: gi|209571234 cell wall protein V [[Clostridium] difficile]

Peptides of selected protein

Standardised Expression Profiles

Section Complete

This leaves the peptides with the same expression pattern assigned with the appropriate protein thus resolving the conflicts.

LC-MS Tutorial - Progenesis LC-MS

Experiment design: AC

Proteins: No filter applied

Accession	Peptides	Conflict	Score	Tags
gi 384359782 (+8)	21	0	1.88E+03	
gi 209571234 (+3)	12	0	1.25E+03	
gi 126699078 (+2)	3 (1)	3	307	
gi 126699128	3 (1)	3	328	
gi 254976383 (+1)	5	0	428	
gi 255656774 (+2)	1	0	136	

Peptides of gi|209571234

#	Score	Hits	Mass	Mass error (p...)	RT (mins)	Charge	Tags	Abundance	Conflict	Peptide
2183	108	4	1821.024	-0.673	38.2	2		9E+04	0	K.V.
1215	64.2	5	2171.126	0.221	59.1	3		1.48E+05	0	D.G.
1508	67	5	1609.903	-0.303	61.1	3		5.68E+04	0	Q.D.
2009	61.3	5	1708.867	0.115	46.2	3		1.38E+05	0	Y.V.
459	116	5	1708.871	2.47	46.1	2		1.54E+06	0	Y.V.
1405	68.6	5	2039.964	-0.31	30.1	3		4.93E+05	0	A.G.

Protein: gi|209571234 cell wall protein V [[Clostridium] difficile]


Peptides of selected protein

Standardised Expression Profiles

Section Complete

**Note:** the number of conflicts you have to resolve will depend on the scope and stringency of the filters you apply at the **Refine Identifications** stage.

**Note:** the default **Protein options** for protein grouping and Protein quantitation are set as shown



With protein grouping switched on protein groups and the additional members are indicated by a bracketed number located after the Accession number. Taking **flagellin subunit** as an example, when the cursor is held over the accession number the group members (9 in total) appear in a tool tip.

Experiment design: AC

Protein: No filter applied Create...

Peptides of gi|260682017: No filter applied Create...

Accession	Peptides	Conflict	Score	#	Score	Hits	Mass	Mass error (p...)	RT (mins)	Charge	Tags	Abundance	Co
gi 254976385 (+1)	11	0	1.06E+03	76	63.9	5	1669.888	-0.0985	42.4	3		4.82E+06	
gi 260682017 (+8)	12	0	1.35E+03	138	103	5	1669.887	-0.663	42.4	2		4.32E+06	
gi 255656776 (+1)								0.407	22.9	2		3.11E+06	
gi 126700407								0.167	38.9	2		5.6E+06	
gi 126697969 (+1)								0.157	38.9	3		3.12E+06	
gi 126697970								0.445	30.5	2		1.75E+06	
gi 126700790 (+1)								-1.19	34.8	3		7.66E+05	
gi 126698640								-1.21	34.8	2		7.54E+05	

↑ Protein: gi|260682017 flagellin subunit [Clostridium difficile CD196]

↓ No protein selected

Peptide Views Protein Resolution

Conflicting proteins

Accession	Peptides	Conflict	Protein Score

Peptides of conflicting protein

#	Score	Hits	Mass	Mass error (p...)	RT (mins)	Charge	Tags	Abundance	Conflict

Protein options... Section Complete

Having performed the conflict resolution with **Group similar proteins** and **Quantitate from non-conflicting features** now switch off the protein grouping.

As grouping is switched off the grouped proteins appear with conflicts to the other group members and the number of unique peptides that are used for quantitation appear in brackets after the peptide number.

**Note:** flagellin has **no unique** peptides (brackets after the peptides field) as they are all present in flagellin subunit protein hence the reason for grouping. As a result all the conflicts are internal to the group.

Experiment design: AC

Protein: No filter applied Create...

Accession	Peptides	Conflict	Score	Tags
gi 5668937	12 (0)	75	1.35E+03	
gi 260682017	12 (0)	75	1.35E+03	
gi 126697810	9 (0)	65	1.07E+03	
gi 10281485	8 (0)	60	867	
gi 357548088	7 (0)	56	821	
gi 10281487	8 (0)	55	887	
gi 73745726	5 (0)	46	662	
gi 209570719	10 (0)	37	816	

Peptides of gi|260682017 No filter applied Create...

#	Score	Hits	Mass	Mass error (p...)	RT (mins)	Charge	Tags	Abundance	Conf
3358	84.7	4	1423.649	-0.431	22.7	2		7.25E+04	
449	93.3	5	1676.838	-1.21	34.8	2		7.54E+05	
1845	99	4	1407.656	0.593	38.7	2		1.28E+05	
477	47.4	5	1692.835	-0.238	20.7	3		5.61E+05	
888	109	5	1692.833	-1.02	20.7	2		4.67E+05	
300	43.8	5	1676.838	-1.19	34.8	3		7.66E+05	
157	125	5	2317.115	0.167	38.9	2		5.6E+06	
175	60.9	5	2317.115	0.157	38.9	3		3.12E+06	

↑ Protein: gi|260682017 flagellin subunit [Clostridium difficile CD196]  
↓ Protein: gi|5668937 flagellin [[Clostridium] difficile]

Peptide Views Protein Resolution

Conflicting proteins for feature 3358

Accession	Peptides	Conflict	Protein Score
gi 5668937	12 (0)	75	1.35E+03
gi 260682017	12 (0)	75	1.35E+03
gi 126697810	9 (0)	65	1.07E+03
gi 10281487	8 (0)	55	887
gi 10281485	8 (0)	60	867
gi 357548088	7 (0)	56	821
gi 73745726	5 (0)	46	662

Peptides of gi|5668937

#	Score	Hits	Mass	Mass error (p...)	RT (mins)	Charge	Tags	Abundance	Conf
3358	84.7	4	1423.649	-0.431	22.7	2		7.25E+04	6
449	93.3	5	1676.838	-1.21	34.8	2		7.54E+05	6
1845	99	4	1407.656	0.593	38.7	2		1.28E+05	6
477	47.4	5	1692.835	-0.238	20.7	3		5.61E+05	6
888	109	5	1692.833	-1.02	20.7	2		4.67E+05	6
300	43.8	5	1676.838	-1.19	34.8	3		7.66E+05	6
157	125	5	2317.115	0.167	38.9	2		5.6E+06	5
175	60.9	5	2317.115	0.157	38.9	3		3.12E+06	5
2288	57.7	3	1386.71	-0.0289	17.8	3		4.61E+04	5

Protein options... Section Complete

Now set the Protein Options back to **Group similar proteins**

Now move to the Review Proteins section by clicking on **Review Proteins** icon on the workflow at the top of the screen.

## Stage 12: Review Proteins

The **Review Proteins** stage opens displaying details for all the proteins. You can now create tags at the protein level by right clicking on the table and selecting **Quick Tags**...

Accession	Peptide count	Confidence score	Anova (p)	Tag	Max fold change	Highest Mean	Lowest Mean	Description
254976387 (+6)	33	3.75E+03	5.33E-06		8.06	C	A	cell surface protein (S-layer precursor protein) [Clostridium difficile QCD-630]
126701103	1	51.7	7.17E-06		6	A	C	ribose-5-phosphate isomerase 2 [Clostridium difficile 630]
54781345 (+1)	5	398	1E-05		8.26	A	C	(R)-2-hydroxyisocaproate dehydrogenase [Clostridium difficile ATCC 9689]
126699971	5	334	1.2E-05		11.8	A	C	<b>thioredoxin 2 (Trx2) [Clostridium difficile 630]</b>
126697690	5	498	3.43E-05		6.01	A	C	ferredoxin/ferredoxin oxidoreductase subunit gamma [Clostridium difficile 630]
126700634	3	48.9	4.38E-05		6.41	A	C	PTS system mannose-specific transporter subunit IIB [Clostridium difficile 630]
126699140	1	158	5.19E-05		12.4	A	C	ferredoxin-NADP(+) reductase subunit alpha [Clostridium difficile 630]
126699940	1	79.3	5.3E-05		4	A	C	transketolase, N-terminal (Sedoheptulose-7-phosphate:D-glyceraldehyde-3-phosphate transferase) [Clostridium difficile 630]
260682017 (+8)	12	1.88E+03	7.52E-05		3.43	A	C	flagellin subunit [Clostridium difficile CD196]
126697752	5	498	3.43E-05			A	C	NAD-specific glutamate dehydrogenase [Clostridium difficile 630]
126698435	1	48.9	4.38E-05			C	A	ABC transporter oligopeptide-family extracellular solute-binding protein [Clostridium difficile 630]
126697684	3	158	5.19E-05			A	C	phosphate butyryltransferase [Clostridium difficile 630]
126697583	1	79.3	5.3E-05		4	A	C	DNA binding protein [Clostridium difficile 630]
126701233	1	44.9	6.3E-05		3.63	A	C	ferredoxin [Clostridium difficile 630]
384359782 (+8)	21	1.88E+03	7.52E-05		3.43	C	A	hemagglutinin/adhesin [Clostridium difficile BI1]

**Selected protein: thioredoxin 2 (Trx2) [Clostridium difficile 630]**  
[View peptide measurements](#)

ArcSinh Normalised Abun...  
 16.0  
15.5  
15.0  
14.5  
14.0  
13.5


Total number of displayed proteins: 135

Section Complete

As an example let us explore Thioredoxin 2.


The table indicates that this protein is most highly expressed in Condition A by 11.8 fold over the lowest condition (C).

To view the corresponding peptide measurements for the current protein either double click on the protein in the table or use the **View peptide measurements** below the table.




**Note:** by selecting all the peptides you can compare the pattern of expression across all the samples allowing you to identify 'atypical' behaviour of peptides assigned to the current protein.

Modified proteins can be located by specifically searching for proteins containing modified peptides. Use the **Back** button to return to the Proteins List and right click on it and select **Modification** from the list of **Quick Tags**.




To find those proteins containing peptides with Oxidation on Methionine residues enter **Oxidation (M)**. This will automatically provide a named tag when you click **Create tag**.





To reduce the table to displaying only these proteins with modified peptides (Oxidation on methionine) use the tag filter to focus on these proteins.



The proteins table will now only display those proteins containing modified peptides.

Accession	Peptide count	Confidence score	Anova (p)	Tag	Max fold change	Highest Mean	Lowest Mean	Description
gij254976387 (+6)	33	3.75E+03	5.33E-06		8.06	C	A	cell surface protein (S-layer precursor protein) [Clostridium difficile QCD-66C]
gij126699971	5	334	1.2E-05		11.8	A	C	thioredoxin 2 (Trx2) [Clostridium difficile 630]
gij126697690	5	491	1.32E-05		6.01	A	C	ferredoxin/ferredoxin oxidoreductase subunit gamma [Clostridium difficile 630]
gij260682017 (+8)	12	1.35E+03	3.33E-05		3.87	A	C	flagellin subunit [Clostridium difficile CD196]
gij384359782 (+8)	21	1.88E+03	7.52E-05		3.43	C	A	hemagglutinin/adhesin [Clostridium difficile B11]
gij126700790 (+1)	8	688	9.71E-05		2.41	C	A	enolase (2-phosphoglycerate dehydratase) (2-phospho-D-glycerate hydro-ly
gij126700129	3	267	0.000199		2.4	A	C	translation inhibitor endoribonuclease [Clostridium difficile 630]
gij126697654	3	215	0.000219		2.53	A	C	30S ribosomal protein S8 [Clostridium difficile 630]
gij54781347	4	341	0.00222		2.47	A	C	2-hydroxyisocaproate-CoA transferase [Clostridium difficile ATCC 9689]

Selected protein: thioredoxin 2 (Trx2) [Clostridium difficile 630]  
[View peptide measurements](#)

Total number of displayed proteins: 9


Section Complete

**Note:** the **Sequence Quick** tag can be used to locate Proteins containing peptides with specific motifs.


Now move to the Protein Statistics section by clicking on **Protein Statistics** icon on the workflow at the top of the screen.

## Stage 13: Protein Statistics

Protein Statistics opens with a Principal Components Analysis (PCA) for all the proteins displayed.




The Multivariate Stats can now be applied to all or subsets of proteins as determined by the current Tag filters (in this example none). Allowing you to identify similar patterns of expression using the Correlation Analysis.



Now move to the **Report** section to report on Proteins and /or peptides.


## Stage 14: Reporting

The **Report Design** stage allows you to select what views you want to include in a report based on the list of **currently selected proteins**.



As an example we will create a report for **only** the proteins containing peptides with Oxidation of Methionine residues enter **Oxidation (M)**.

1. First reduce the proteins to report on by selecting the '**Modification with Oxidation (M)**' tag. In this example it reduces the number of proteins in the table to 9.
2. Expand the various Report Design options (by default they are all selected)
3. Un-tick as shown below
4. Click **Create Report**



This opens a dialog to allow you to save the report, after which it will be opened in the form of a web page.

Click on the **Accession No.** in the proteins section of the Report and this will take you to the Assigned peptides for this protein

LC-MS Tutorial

Experiment: LC-MS Tutorial  
Report created: 18/02/2013 09:39:30

Proteins

Protein building options  
Protein grouping **Group similar proteins**  
Protein quantitation **Using only features with no protein conflicts**

Accession	Peptides	Score	Anova (p)	Fold	Tags	Description	Average Normalised Abundances	
							A	C
gi 254976387	33	3752.13	5.33e-006	8.06		cell surface protein (S-layer precursor protein) [Clostridium difficile]	5.35e+007	4.32e+008
gi 384359782	21	1879.25						
gi 260682017	12	1352.16						
gi 126700790	8	688.34						
gi 126697690	5	490.99						
gi 54781347	4	341.12						
gi 126699971	5	333.68						
gi 126700129	3	266.82						
gi 126697654	3	214.79						

gi|126699971

thioredoxin 2 (Trx2) [Clostridium difficile 630]  
5 peptides

Sequence	Feature	Score	Hits	Mass	Charge	Tags	Conflicts	Modifications	In quantitation	Average Normalised Abundances	
										A	C
DDATVPNIENMIK	1429	60.28	5	1458.7019	2		0		yes	2.29e+005	1.37e+004
VDEVTKDDATVPNIENMIK	1157	66.32	4	2130.0505	2		0		yes	5.70e+005	9.09e+004
VDEVTKDDATVPNIENMIK	366	45.25	4	2130.0508	3		0		yes	1.40e+006	1.38e+005
VDEVTKDDATVPNIENMIK	1730	44.06	2	2146.0465	3		0	[17] Oxidation (M)	yes	2.17e+005	1.74e+004
VLGLPTMAIYK	263	69.09	5	1204.6890	2		0		yes	2.29e+006	1.40e+005
VLGLPTMAIYK											

Accession gi|126699971

Description thioredoxin 2 (Trx2) [Clostridium difficile 630]  
Peptides 5  
Score 333.68  
Anova 1.20e-005  
Fold 11.77

- Anova p-value ≤ 0.05
- Max fold change ≥ 2
- Modification with Oxidation (M)

Tags

- Most abundant
- Anova p-value ≤ 0.05
- Max fold change ≥ 2
- Significantly up in C
- Up regulated in C
- Up regulated in A

Having closed the report it can be reopened by double clicking on the saved html file.



**Note:** you can also copy and paste all or selected sections of the report to Excel and/or Word.

## Exporting Protein Data



Protein data can be exported in a csv file format. You can either export the **Protein and/or peptide measurements** using the options in the File Menu which are available at the Review Proteins stage.

As an example of Data export use the Tag filtered set from the previous section. Where you are only going to export measurements for those proteins that have Oxidised Methionine residues.

First set the tag filter as shown below. Then select **Export Protein Measurements** from the **File** menu.



The Export Protein Measurements dialog opens. Select the required fields and click OK. **Save** the file and then open the exported data file using the dialog that opens.




Excel will open displaying the exported protein measurements

	A	B	C	D	E	F	G
1						Normalized abundance	
2						A	A2
3	Accession	Peptide count	Anova (p)	Max fold change	Description	A1	A3
4	gi 254976387;gi 112	33	5.33E-06	8.063180835	cell surface protein (S-layer precursor protein) [Clostridium	55837799	57754081
5	gi 126699971	5	1.20E-05	11.7737882	thioredoxin 2 (Trx2) [Clostridium difficile 630]	4957161	5277038
6	gi 126697690	5	1.32E-05	6.011328857	ferredoxin/ flavodoxin oxidoreductase subunit gamma [Clos	1778402	2084321
7	gi 260682017;gi 102	12	3.33E-05	3.872261726	flagellin subunit [Clostridium difficile CD196]	27038907	24794152
8	gi 384359782;gi 209	21	7.52E-05	3.430214543	hemagglutinin/adhesin [Clostridium difficile BI1]	2614529	2673778
9	gi 126700790;gi 296	8	9.71E-05	2.414242661	enolase (2-phosphoglycerate dehydratase) (2-phospho-D-g	404610.9	387979.1
10	gi 126700129	3	0.000198776	2.398837692	translation inhibitor endoribonuclease [Clostridium difficile	5557391	5446163
11	gi 126697654	3	0.0002192	2.532221342	30S ribosomal protein S8 [Clostridium difficile 630]	1014579	1149285
12	gi 54781347	4	0.002217665	2.469813207	2-hydroxyisocaproate-CoA transferase [Clostridium difficile	1581089	1205546

**Note:** where there are multiple group members the other accession numbers are also exported.

At each stage in the Work flow there are a number of Export and Import options available from the **File** Menu. This includes the option to **Import Additional Protein Data** which can be used to increase the Protein meta data and also be used to sort the existing tabular data.



## Creating an Inclusion list


Inclusion lists can be used to try to increase the number of identified peptides you have in your experiment. They are used to control your mass spectrometer to try and concentrate the collection of MS/MS scans from specified m/z – retention time positions. Runs made using the inclusion list can then be imported into the Progenesis LC-MS and the extra MS/MS scans added to the experiment.

As an example we are going to create an inclusion list for all the features that show a Significant difference between groups A and C (Anova p<0.05) and have **no MS/MS spectra**.




First return to **Review Peak Picking** using the Workflow icons.


With no filters applied right click on a feature in the table, select **Quick tags** and click on **No MS/MS data**.




Filter the table so that it is only showing features with an **Anova p-value ≤ 0.05** and **No MS/MS data**. To do this click **Create** and drag the two tags on to the Show features that have all of these tags and click OK.



Select all of the features that are displayed, right click and create a New tag called **Inclusion 1**




Call the new tag Inclusion\_1.








Now use the new tag to filter the table to display only those features that show a Significant Change and **DO NOT** have any MS/MS spectra.



Then select **Export Inclusion List...** from the file menu and select the appropriate format.



Finally save the file to an appropriate location



**Note:** with certain MS machines it is possible to widen the retention time windows being used, this can be controlled using the following dialog.

If you require further information on the inclusion list file formats then click **Help**.

**Note:** The new LC-MS runs can then be added to the original experiment to increase the MS/MS coverage using the **Add files** facility at the Data Import Stage.

## Congratulations!

This document has taken you through a complete analysis using Progenesis LC-MS, from Alignment through Analysis to generating lists of interesting features using powerful Multivariate Statistical analysis of the data.

Hopefully our example has shown you how this unique technology can deliver significant benefits with

- Speed
- Objectivity
- Statistical Power

If you would like to see the benefits of running Progenesis SameSpots using your own runs and explore the Progenesis LC-MS workflow please go to Appendix 3: Licensing Runs (page 60).


## Appendix 1: Stage 1 Data Import and QC review of LC-MS data set

You can use your own data files, either by directly loading the raw files (Waters, Thermo, Bruker and Agilent) or, for other Vendors, convert them to mzXML or mzML format first.

To create a new experiment with your files select **New** give your experiment a name. Then select data type, the default is 'Profile data'.


**Note:** if you have converted or captured the data as centroided then select Centroided data and enter the Resolution for the MS machine used.

Click **Create experiment** to open the LC-MS Data Import stage of the workflow.



Select the 'Import Data file format', in this example they are mzXML files

Then locate your data files using **Import...**




Locate and select all the Data files (A1 to C3).

On loading the selected runs your data set will be automatically examined and the size of each file will be reduced by a 'data reduction routine', which reduces the data by several orders of magnitude but still retains all the relevant quantitation and positional information.


**Note:** For a large number of files this may take some time.

**Note:** as the loading process starts you can also start the automatic alignment before the loading has completed. This is a 2 stage process that involves the selection of an Alignment Reference (either automatically or manually then the automatic alignment of all your runs to this Reference run.

Click **Start alignment process** to start the automatic alignment of your runs.



Progenesis LC-MS provides three methods for choosing the alignment reference run, as seen below:



### 1. Assess all runs in the experiment for suitability


This method compares every run in your experiment to every other run for similarity.

The run with the greatest similarity to all other runs is chosen as the alignment reference.

If you have no prior knowledge about which of your runs would make a good reference, then this choice will normally produce a good alignment reference for you. This method can take a long time

### 2. Use the most suitable run from candidates that I select

This method asks you to choose a selection of reference candidates, and the automatic algorithm chooses the best reference from these runs.



When you have some prior knowledge of your runs suitability as references:

runs from pooled samples

runs for one of your experimental conditions will contain the largest set of common peptides.


### 3. Use this run

This method allows you to manually choose the reference run.

Manual selection gives you full control, but there are a couple of risks to note:


- If you choose a pending run which subsequently fails to load, alignment will not be performed.
- If you choose a run before it fully loads, and it turns out to have chromatography issues, alignment will be negatively affected (for this reason we recommend that you let your reference run fully load and assess its chromatography before loading further runs).

Once you have selected how to handle the choice of Reference run you will now be asked if you want to Align your runs automatically or manually.



Select automatically and click finish.

The Alignment process starts with the automatic selection of C1 as the reference




Once the Reference run has been chosen the automatic alignment is performed. As the whole process proceeds you get information on what stage has been performed and also the % of the process that has been completed.

**Note:** At this stage you have the option to Review the Chromatography or go straight to the review of the Automatic Alignment of your data.

### Review Chromatography

Each data file appears as a 2D representation of the run. At this stage you will be warned if any of the data files have been 'centroided' during the data acquisition and conversion process.

You can delete run(s) by left clicking on the run in the list.



At the Import Data stage you can examine the quality of the imported runs using the 2D representation of the runs

**Import Data**

- 1 Import your run data**  
Select one of the available data formats then click the Import button:  
Format:    
[About this data format](#) [Download others](#)
- 2 Start the alignment process**  
While your runs are importing, click the button below to:
  - Select [alignment reference](#) candidates ✓
  - Determine the best of the candidates ✓
  - Align all runs to that reference run ✓
- 3 Review the chromatography**  
Look at all of the runs in the list below, checking for any [sample-running problems](#) that might affect analysis. Right-click to remove any runs that have significant problems.

Imported runs:

A1 A2 A3  
C1 C2 C3

**C1**

**About this run**

- MS/MS count: 9,542
- Peak count: 1,182,000
- Total ion intensity: 3.643e+009
- Masked areas : none

**Alignment reference**  
This run is being used as the experiment's alignment reference.  
If you want to use a different run as the alignment reference, you'll need to discard any analysis and restart the alignment process:

Zoom:

**Note:** details of the current run appear on the top right of the view.

Once you have reviewed the imported runs click on Review Alignment or Section Complete to move forward to the Review Alignment Stage.

**Note:** you will be offered the automatic alignment if you have not performed it automatically already.

Now move to the next stage in the workflow (page 9 in this user guide) by clicking **Section Complete**.





## Appendix 2: Stage 1 Data QC review and addition of exclusion areas

During the process of Data QC you may identify areas of the raw data for a particular run that appear 'noisy' yet still have identifiable 'isotopic patterns'.

For example if the run is part of a 'replicate set' of runs it is possible to exclude such areas on the noisy run by applying a mask to the area. By doing so this area is excluded during the initial part of the detection process in order that it does not 'interfere' with the detection of the features in the replicate group.

To do this select Masked areas from selected run on the bottom left of the screen.

Drag out an area over the noisy part of the run to create the mask.




**Note:** Click **Done** to return to the **Import Data** view where you can zoom into the masked where you will see the isotopic features in the noise.

**Note:** if the level of noise is high and affecting many of your runs a preferred approach would be to re-optimize the chromatography to improve the levels of noise in your data.



## Appendix 3: Licensing runs (Stage 3)

When setting up a **New experiment** if you are evaluating Progenesis LC-MS with unlicensed runs then the licensing page will open after **Import Data section**.



**If you already have a programmed dongle attached to your machine then the License Runs page will not appear.**


To use this page to License your Runs you must first either obtain an **'Evaluation' Licence Code** from a **Nonlinear Sales Person** or **purchase a licence code directly from Nonlinear**.

Each code will allow you to license a set number of runs.


The runs in your experiment will be listed as shown below.

To activate license(s) for the selected runs enter the code in the space provided and click **Use Licence code**.

**Note:** you will need an internet connection to use this method.



A message confirming successful installation of your licences will appear.



Click **OK**, the view will update and Alignment, the next stage in the workflow, will open with the licensed files.

## Appendix 4: Manual assistance of Alignment

### Approach to alignment

To place manual alignment vectors on a run (C2 in this example):


1. Click on Run C2 in the **Runs** panel, this will be highlighted in green and the reference run (C1) will be highlighted in magenta.
2. You will need approximately 5 - 10 **alignment vectors** evenly distributed from top to bottom of the whole run.
3. First either drag out a focus area or set the focus area to **8** or **16** using the Focus grid size on the bottom left of the screen.




Click on an area (see below) in the **Whole Run** view (C) to refocus all the windows. Adjust Contrast as required.

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/>	0	79.0%
A2	<input checked="" type="checkbox"/>	0	43.7%
A3	<input checked="" type="checkbox"/>	0	37.7%
C1	<input checked="" type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/>	0	38.6%
C3	<input checked="" type="checkbox"/>	0	82.3%

- Click and hold on a green feature in Window A as shown below.




- As you are holding down the left mouse button drag the green feature over the corresponding magenta feature of the reference run. The red circle will appear as shown below indicating that a positional lock has been found for the overlapping features.



**Note:** as you hold down the mouse button, window B zooms in to help with the alignment.


- On releasing the left mouse button the view will 'bounce' back and a red vector, starting in the green feature and finishing in the magenta feature will appear.






**Note:** an incorrectly placed vector is removed by right clicking on it in the **Vector Editing** window

7. Now click **Show Aligned** on the top tool bar to see the effect of adding a single vector.



8. Adding an additional vector will improve the alignment further. **Note** this time as you click to add the vector it 'jumps' automatically to the correct position using the information from the existing alignment vector.



Repeat this process moving the focus from top to bottom on the **Whole Run** view

**Review Alignment**  
Align peptide ions to compensate for drifts in retention time by dragging them up or down in the Vector Editing window.

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/> <input type="checkbox"/>	0	79.0%
A2	<input checked="" type="checkbox"/> <input type="checkbox"/>	0	43.7%
A3	<input checked="" type="checkbox"/> <input type="checkbox"/>	0	37.7%
C1	<input checked="" type="checkbox"/> <input type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/> <input type="checkbox"/>	5	96.4%
C3	<input checked="" type="checkbox"/> <input type="checkbox"/>	0	82.3%

Contrast:

Ion maps:  Alignment target  Run being aligned

Focus grid size:  Make box square

Alignment quality: ■ Good ■ OK ■ Needs review

Section Complete


**Note:** the number of vectors you add is recorded in the **Runs** table also with each vector addition the Score and alignment quality updates. This can help guide the number of manual vectors you need to add before applying the automatic alignment.

10. Repeat this process for all the runs to be aligned.


The number of manual vectors that you add at this stage is dependent on the misalignment between the current run and the Reference run. In many cases only using the Automatic vector wizard will achieve the alignment.

Also the 'ease' of addition of vectors is dependent on the actual differences between the LC-MS runs being aligned.





11. Then select Automatic Alignment to bring up the Automatic Alignment dialog and click **OK**. The automatic alignment process will begin, using the manual vectors you have added to aid in automatic vector placement.



**Note:** the tick boxes next to the 'Run name control' which control whether vectors will be generated for each run.

To review the vectors, automatic and manual return to page 11

## Appendix 5: Within-subject Design

To create a **Within-subject Design** for your data set select this option on the **Experiment Design Setup** page and enter the name of the design.

In this example there are 3 Subjects (i.e. patients A, B and C) who have been individually sampled: Before(1), During (2) and After (3) treatment

Which experiment design type do you want to use for this experiment?

Between-subject Design

**Within-subject Design**

Do samples from a given subject appear in only one condition? Then use the between-subject design.

To set up this design, you simply group the runs according to the condition (factor level) of the samples. The ANOVA calculation assumes that the conditions are independent and therefore gives a statistical test of whether the means of the conditions are all equal.

Have you taken samples from a given subject under different conditions? Then use the within-subject design.

Note: you must have a sample from every subject for every condition to


assumption of independence. With a repeated measures ANOVA individual differences can be eliminated or reduced as a source of between condition differences (which helps to create a more powerful test).

The within-subject design can be thought of as an extension of the paired-samples t-test to include comparison between more than two repeated measures.

	Before	During	After
Patient X	X1	X2	X3
Patient Y	Y1	Y2	Y3
Patient Z	Z1	Z2	Z3

When the design page opens use the **Add Subject** and **Add Condition** buttons to create the matrix that fits your experimental design, over typing the names as required.

Then Drag and drop the Samples on to the correct 'cell' of the matrix.



You can create additional Experimental Designs using the New tab

All of these Experimental Designs are available at all the following stages in the LC-MS workflow.

## Appendix 6: Power Analysis (Peptide Stats)

Power analysis is a statistical technique that is used to gauge how many replicates are needed to reliably see expression differences in your data. It is available through the Progenesis Stats section of the workflow.

To perform a power analysis of the data click on **Ask another question** at the top of the table in the Progenesis Stats screen. A selection of 3 tools will appear in the form of questions.

**Principal Components Analysis**  
Are there any outliers in my data?  
Does my data cluster according to my experimental conditions?

**Correlation Analysis**  
Group my features according to how similar their expression profiles are.

**Power Analysis**  
How many replicates should I run?  
What is the power of my experiment?


Select the option

*'How many replicates should I run and what is the power of my experiment?'*

It answers this question by informing you:

*'How many replicates you need so that at least 80% of your features with a power >0.8'*

Using the **Significant p<0.05 features (7965)**, as an example, view the power analysis.



This is displayed graphically showing that 78.6% of the 7965 features have a power of 80% or that 4 replicates would give you 96.4% of your data with power > 0.8.

- **The power of a statistical test** reflects our confidence in the experimental data's ability to find the differences that do actually exist
- The power is expressed as a percentage, where 80% power is an accepted level, therefore allowing you to assess the number of sample replicates that would be required to achieve a power of 80%.

## Appendix 7 Search engine parameters (Stage 9) Mascot

The parameters applied to the Mascot search that yielded the search results used in this user guide are shown below:

MASCOT MS/MS Ions Search			
<b>Your name</b>	andy.borthwick	<b>Email</b>	andy.borthwick@nonlinear.com
<b>Search title</b>	Search_full_LC-MS_4.1		
<b>Database(s)</b>	ApoE MSDB NCBIInr NIST_Yeast_Sigma SwissProt	<b>Enzyme</b>	Trypsin
		<b>Allow up to</b>	1 missed cleavages
		<b>Quantitation</b>	None
<b>Taxonomy</b>	... Firmicutes (gram-positive bacteria)		
<b>Fixed modifications</b>	Carbamidomethyl (C)	> <	mTRAQ:13C(6)15N(2) (N-term) mTRAQ:13C(6)15N(2) (Y) NIPCAM (C) Oxidation (HW) Propionamide (C) Pyridylethyl (C) Pyro-carbamidomethyl (N-term C) Sulfo (STY) TMT2plex (K) TMT2plex (N-term) TMT6plex (K)
	Display all modifications <input type="checkbox"/>		
<b>Variable modifications</b>	Oxidation (M) Phospho (ST) Phospho (Y)	> <	
<b>Peptide tol. ±</b>	9 ppm # <sup>13</sup> C 0	<b>MS/MS tol. ±</b>	0.6 Da
<b>Peptide charge</b>	2+	<b>Monoisotopic</b>	<input checked="" type="radio"/> Average <input type="radio"/>
<b>Data file</b>	D:\Customer Data\LC-MS\LC-ME Browse...		
<b>Data format</b>	Mascot generic	<b>Precursor</b>	m/z
<b>Instrument</b>	ESI-TRAP	<b>Error tolerant</b>	<input type="checkbox"/>
<b>Decoy</b>	<input type="checkbox"/>	<b>Report top</b>	AUTO hits
<b>Start Search ...</b>		<b>Reset Form</b>	

Database : NCBIInr (circa 02/13) was used with the Taxonomy restriction set to Firmicutes

Fixed modifications: Carbamylation(C) and variable modifications: Oxidation (M), Phospho (ST) and Phospho (Y)

Peptide Tol: 9ppm

Instrument: ESI-Trap